



Pressure Care Made Eezi



RECLINERS

- *The latest and easiest solution to eliminate the risk of Pressure Sores.*
- *Wicks away moisture.*
- *Breathable.*
- *Totally eliminates shearing forces by minimising the friction between the skin and pad.*
- *Can be washed and bleached.*
- *Bed Sheet Size is 100cm x 80cm.*
- *Weight Capacity:
2 Layer - 12 Stone
3 Layer - 17 Stone
4 Layer - 25 Stone*

DESIGN * CRAFT * QUALITY

Instructions

Wash the Treat-Eezi pad before use. Place directly under the pressure sore or an area at risk of developing a pressure sore. Best results occur when the sore is directly in contact with the product. If a garment or undergarment is worn, ensure it is thin as possible and allows air circulation (e.g a thin cotton shirt or thin cotton undergarment.) Do not place any other pad, bedding or cushions over the Treat-Eezi pad. Skin creams can be used in conjunction with the Treat-Eezi pad.

Wash using a pH neutral detergent. Decontaminate using Sodium Hypochlorite solution of 1000ppm in line with HSG9518 and BHTA guidelines. Ensure product is completely dry before patient use. Dry naturally or on a cool tumble setting.

DO NOT IRON & DO NOT USE FABRIC SOFTENERS

Disclaimer - considering individual circumstances may be appropriate for the use of the products standard medical care must be followed at all times. These products are not intended to replace and procedure or guidelines in place for the prevention of pressure wounds.

Warranty - this product is guaranteed to be free from defects in materials and workmanship under normal use for a period of one year.

To find out more about our Pressure Care products and various Seating Ranges contact us direct or visit us online.

Recliners, Old Griffin Field, Windsor Street, RCT, South Wales, UK, CF41 7JJ
info@reclinersltd.co.uk | 01443 420851 / 07977494774

www.reclinersltd.co.uk

